

# MAXIM<sup>®</sup>

## Lemon Detergent Disinfectant

### Features:

**Disinfectant & Fungicide**  
**Animal Premise Virucidal**  
**Mildewstat, Virucide**

### Description

This product is a no-rinse disinfectant cleaner that disinfects, cleans and deodorizes in one labor-saving step. This product has also been designed as an effective disinfectant for veterinarian and animal care facilities. This product can be used on washable hard, non-porous surfaces such as: glass surfaces, metal, plastic, sealed stone, and vinyl.

### Directions

#### FOR USE AS A ONE-STEP, GENERAL, HOSPITAL, MEDICAL DISINFECTANT, FUNGICIDE, VIRUCIDE\*, CLEANER:

- 1 Pre-clean heavily soiled areas.
- 2 Apply use solution of 2 oz. of this product per gal. of water to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, mechanical spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.
- 3 Treated surfaces must remain wet for 10 minutes.
- 4 Allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.
- 5 Prepare a fresh solution daily or when visibly dirty.



### Applications

#### This product is for use in:

- > Hospitals and healthcare facilities
- > Veterinary clinics and animal premises
- > Schools

### Specifications

**Color** ..... Yellow  
**Fragrance** ..... Lemon  
**pH Range** ..... 11.5 - 12.0  
**Dilution** ..... 2 oz./gal.  
**Density** ..... 8.59

### Packaging

4/1 gallon	047000-41
5 gallon	047000-05

### Product Registration No.

EPA: 10324-93-45745

ESL 012615