

# MAXIM<sup>®</sup>

## Maxim Hi-Q

### Features:

*Hospital-use disinfectant*

*Deodorizes*

*Concentrated*

### Description

This product can be used on hard, non-porous surfaces such as: floors, walls, tables and chairs. Kills Pandemic 2009 H1N1 Influenza A Virus, (formerly called swine flu).

### Directions

#### FOR USE AS A GENERAL, HOSPITAL, MEDICAL DISINFECTANT, VIRUCIDE\*, CLEANER:

- 1 Pre-clean heavily soiled areas.
- 2 Apply use solution of 0.5 oz. of this product per gal. of water to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, auto scrubber, or mechanical spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.
- 3 Treated surfaces must remain wet for 10 minutes.
- 4 Allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.
- 5 Prepare a fresh solution daily or when visibly dirty.



### Applications

#### This product is for use in:

- > Hospitals
- > Public restrooms
- > Schools

### Specifications

Color .....	Blue
Fragrance .....	Floral
pH Range .....	6.0 - 6.5
Dilution .....	1/2 oz./gal.
Density .....	8.42

### Packaging

4/1 gallon	047300-41
30/2 oz. packet	047300-60

### Product Certifications

EPA: 10324-108-45745

E15